


# Series 2000XL

**3-sectioned multipurpose table with electric motorized height and inclination adjustment** ... for examination, rectoscopy, gynaecology ...

- Model 2648-00XLE white powder coated (RAL 9010)
- Model 2648-03XLE white aluminium powder coated (RAL 9006)
- Model 2648-04XLE grey aluminium powder coated (RAL 9007)



Electric motorized  
height adjustment



Model 2648-00XLE with options:  
Loss of basic mobility (1182); movable on twin castors with central lock (1180N), paper roll holder (1265K); upper picture additionally:  
Leg supports acc. to Goepel (1520); Upholstery colour: 64-715 anthracite

## Basic equipment:

- Movable on comfort twin castors, Ø 100 mm, each with lock
- Hand grips on both sides of the seat section
- Drawable basin (stainless steel) under seat, basin support: HPL panel
- Including clamps for leg supports (leg supports = available as an option)
- Electromotive drive control of height adjustment by foot switch
- Electromotive drive control of inclination adjustment by hand switch, inclination adjustment to -15°
- Integrated safety function for automatic blocking of electromotive drive (art. no. 1999, mandatory equipment)
- Freewheeling downwards, automatic declutch in case of resistance
- Dewert-hygiene-equipment



## Upholstery:

Lying surface 3-sectioned:

- Head part 30° positive adjustable by 2 metal rachets
- Length head part: 550 mm
- Middle part 50° positive adjustable by rack system
- Length middle part: 930 mm
- Seat section with half-round cutout, manually inclinable (20°)
- Length seat section: 450 mm
- Thickness of upholstery: 80 mm

## Technical data:


- Electronical height adjustment: 650 mm to 1090 mm in 22 sec. by AC-motor with a pressure capacity of 10000 N
- Width of lying surface: 650 mm
- Length of lying surface: 1950 mm
- Load capacity: 220 kg; according to IEC 60601 tested with fourfold static load



## Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



Available as an option (see pricelist for additional options)

## Mobility



Art. no. 1182  
Without basic mobility



Art. no. 1180N  
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm

## Nose opening



Art. no. 1610  
Nose opening, in head part

Art. no. 1620  
Upholstery part to close

## Additional motor control options



Art. no. 1920E  
Hand switch (low voltage), with magnet on the back to fasten



Art. no. 1910E  
Foot switch (low voltage)



Art. no. 1940E  
1 foot switch fixation mounted at underframe



Art. no. 1950E  
Foot switch rail, lateral, on both sides, security adjustment

## Paper roll holder



Art. no. 1265K  
Paper roll holder, mounted under head part (not protruding), for width of lying surface 700 mm

## Head part options



Art. no. 1655  
Head part stepless adjustable by gas spring from 0° to +60°

## Side guards



Art. no. 1040  
1 side guard tiltable, chromium-plated, completely made of steel tube, length: 700 mm



Art. no. 1005  
1 side guard to fit on side rail, with attachment clamp, chromium-plated, length: 800 mm

## Add-on elements



Art. no. 1410  
1 piece of side rail, length: 350 mm



Art. no. 1430  
1 attachment clamp for side rail, suitable for pipe Ø 18 mm



Art. no. 1450  
1 Infusion pole for side rail, without attachment clamp

## Leg supports



Art. no. 1510  
1 pair of leg supports



Art. no. 1520  
1 pair of leg supports according to Goepel

## Upholstery options

Art. no. 1770  
Soft covering for pleasant resting comfort, for width 700 mm, thickness 80 mm

Art. no. 1772  
Upholstery medium-soft (cold foam cover), for width 700 mm, thickness 70 mm

