

Series 2000XL - with continuous frame

2-sectioned versatile table with electric motorized height adjustment

... for examination, treatment, therapy ...

Model 2600-00XLE and 2605-00XLE white powder coated (RAL 9010)

Model 2600-03XLE and 2605-03XLE white aluminium powder coated (RAL 9006)

Model 2600-04XLE and 2605-04XLE grey aluminium powder coated (RAL 9007)



Both pictures:
Model 2600-05XLE
with options: Without standard
mobility (1182), movable
on twin castors with central
lock (1180N), paper roll holder
(1265K);
upholstery colour 64-712 slate
grey


Basic equipment:

- Movable on comfort twin castors Ø 100 mm, each with total lock
- Electromotive drive control by foot switch (option: hand switch)
- Integrated safety function for automatic blocking of electromotive drive
- Freewheeling downwards, automatic declutch in case of resistance
- Dewert-hygiene-equipment

Upholstery:

Lying surface 2-sectioned:

- Head part upwards (+30°) adjustable by 2 metal rachets
- Length head part: 550 mm
- Fixed foot part
- Length foot part: 1400 mm
- Upholstery thickness: 80 mm



Technical data:

- Electric motorized height adjustment: 500 mm to 940 mm in 22 sec. by AC-motor with a pressure capacity of 10000 N
- Width of lying surface: Model 2600-..XLE: 700 mm | Model 2605-..XLE: 800 mm
- Length of lying surface: 1950 mm
- Load capacity: 250 kg; according to IEC 60601 tested with fourfold static load

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology




Available as an option (see pricelist for additional options)

Mobility



Art. no. 1182
Without standard mobility



Art. no. 1180N
Movable on twin castors,
Ø 100 mm with central
lock, basic height +20 mm



Art. no. 1185
Movable on comfort twin
castors, Ø 125 mm with
central lock,
basic height + 50 mm



Art. no. 1045N
Wheel-lifting/fixing
mechanism with 4 foot levers
for controlled lowering

Additional motor control options



Art. no. 1920E
Hand switch (low voltage),
with magnet on the back to fasten



Art. no. 1910E
Foot switch (low voltage)



Art. no. 1940E
1 foot switch fixation mounted at
underframe



Art. no. 1950E
Foot switch rail, lateral, on both sides,
security adjustment



Art. no. 1960E
Foot switch rail, surround,
security adjustment

Paper roll holder



Art. no. 1265K / 1280K
Paper roll holder,
mounted under head part
(not protruding),
for width 650 mm / 800 mm

Head part options



Art. no. 1655
Head part stepless
adjustable by gas spring
(0° to +55°)

Nose opening



Art. no. 1610
Nose opening, in head part

Art. no. 1620
Upholstery part to close

Wall protection bumpers



Art. no. 1170
1 wall protection bumper Ø 150 mm;
in case of side rail 1080: art. no. 1171,
Ø 200 mm

Side guards



Art. no. 1080
1 side guard lateral lowerable, chro-
mium-
plated, completely made of steel tube,
length: 1060 mm



Art. no. 1040
1 side guard tiltable, chromium-plated,
completely made of steel tube,
length: 700 mm

Add-on elements



Art. no. 1410
1 piece of side rail,
length: 350 mm



Art. no. 1430
1 attachment clamp for
side rail, suitable for pipe:
Ø 18 mm



Art. no. 1440
1 arm rest for side rail,
without attachment clamp



Art. no. 1450
Infusion pole, without
attachment clamp