

Heavy-weight tables


2-sectioned heavy-weight table with fixed height

... for examination, treatment, therapy ...

Model 208-00 and 209-00 white powder coated (RAL 9010)

Model 208-03 and 209-03 white aluminium powder coated (RAL 9006)

Model 208-04 and 209-04 grey aluminium powder coated (RAL 9007)



Available as an option (see pricelist for additional options)

Mobility



Art. no. 1189
Movable on heavy-weight twin castors with central lock, Ø 125 mm

Paper roll holder



Art. no. 1210 / 1220
Paper roll holder, for width of lying surface 650 mm / 800 mm

Add-on elements



Art. no. 1810
1 push handle

Side guards



Art. no. 1080
1 side guard foldable, chromium-plated, completely made of steel tube, length: 1060 mm

Head part options



Art. no. 1610
Nose opening

Art. no. 1620
Upholstery part to close

Upholstery options

Art. no. 1770 / 1780
Soft covering for pleasant resting comfort, for width 700 mm / 800 mm, thickness 80 mm



Both pictures::
Model 209-00
with options:
Movable on heavy-weight twin castors with central lock (1189),
paper roll holder (1220);
Upholstery colour: 64-416 black

Basic equipment:

- Table not movable, including height adjustment screws to compensate for uneven ground
- Extremely high load capacity: 400 kg
- Dewert-hygiene-equipment



Upholstery:

Lying surface 2-sectioned:

- Head part stepless adjustable to + 40° by gas spring
- Length head part: 550 mm
- Fixed foot part
- Length foot part: 1400 mm
- Thickness of upholstery: 80 mm

Technical data:

- Height: 650 mm
- Width of lying surface:
Model 208-.. (one longitudinal brace): 650 mm
Model 209-.. (two longitudinal braces): 800 mm
- Length of lying surface: 1950 mm
- Load capacity: 400 kg; according to IEC 60601 tested with fourfold static load



Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology

