

# USEFUL DEVICES FOR DAILY PRACTICE

2023.1

**ZD**  
**DEWERT**

Präzisions-  
Stahlrohrmöbel

MULTI PURPOSE TABLES

PAGE  
2

STOOLS

PAGE  
3

BOLSTER,  
CUSHION, CUBE

PAGE  
6

PAGE  
3

PAGE  
7

STEPS

UPHOLSTERY  
PROTECTION  
AND CARE

German  
manufacture

DURABLE

ROBUST

EASY-CARE



STOOLS



MULTI PURPOSE  
TABLES



STEPS



BOLSTER,  
CUSHION,  
CUBE

# Useful devices for daily practice

## Multi purpose tables



Frame colours  
(left to right:  
-00, -03, -04)

- Model 325N-00 white powder coated (RAL 9010)**
- Model 325N-03 white aluminium powder coated (RAL 9006)**
- Model 325N-04 grey aluminium powder coated (RAL 9007)**

### Equipment:

- Movable on 4 castors Ø 65 mm
- 2 table-tops made of HPL, thickness 8 mm, white
- Big table top (580 x 400 mm) manual tiltable
- Small table top (150 x 200 mm) fixed
- 1 bag-hook under big table top

### Measures:

- Manual height adjustment: 750 - 1.100 mm  
Option: Table lifting unforced by gas spring (art. no. 325-55)
- Measures overall: L: 740 mm x W: 400 mm
- Weight: 11 kg



- Model 327-00 white powder coated (RAL 9010)**
- Model 327-03 white aluminium powder coated (RAL 9006)**
- Model 327-04 grey aluminium powder coated (RAL 9007)**

### Equipment:

- Movable on 4 castors Ø 65 mm
- 1 table-top, 760 x 420 mm, made of HPL, thickness 8 mm, white
- 1 dustbin, white plastics, removable

### Measures:

- Manual height adjustment: 720 - 1.070 mm  
Option: Table lifting unforced by gas spring (art. no. 327-55)
- Measures overall: L: 760 mm x W: 420 mm
- Weight: 13 kg



- Model 329-00 white powder coated (RAL 9010)**
- Model 329-03 white aluminium powder coated (RAL 9006)**
- Model 329-04 grey aluminium powder coated (RAL 9007)**

### Equipment:

- Movable on 4 castors Ø 65 mm, 2 castors with lock
- 2 table-tops, made of HPL, thickness 8 mm, white
- Upper table-top (measures: 650 x 400 mm) height adjustable: 720 mm - 1.020 mm
- Lower table-top (measures: 580 x 400 mm) fixed

### Measures:

- Manual height adjustment of upper table top: 720 mm bis 1.070 mm  
Option: Table lifting unforced by gas spring (art. no. 329-55)
- Measures overall: L: 650 mm x W: 400 mm
- Weight: 15 kg





# Useful devices for daily practice



## Bolster, cushion, cube

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



\*Add colour number



Art. no. 1113-..\*  
Half roll bolster,  
measures (L. x W. x H.):  
570 x 200 x 95 mm



Art. no. 1211-..\*  
Wedge shaped bolster,  
measures (L. x W. x H.):  
400 x 400 x 100 - 10 mm



Art. no. 1112-..\*  
Knee roll,  
measures: L. 570 x Ø 200 mm



Art. no. 1111-..\*  
Neck roll,  
measures: L. 570 x Ø 150 mm



Art. no. 1314-..\*  
Cushion, velcro fixation,  
ticking with polyester filling,  
measures (L. x W.): 400 x 400 mm



Art. no. 1015\*  
Cube, measures (L x W x H):  
400 x 300 x 600 mm

Art. no. 1016\*  
Cube, measures (L x W x H):  
500 x 450 x 400 mm

## Stool with fixed height

**Model 400-00 white powder coated (RAL 9010)**

**Model 400-03 white aluminium powder coated (RAL 9006)**

**Model 400-04 grey aluminium powder coated (RAL 9007)**

- Not movable
- Frame made of square pipe 30 x 30 x 2 mm
- Seat 420 x 420 mm, moulded upholstery, thickness: 65 mm
- Seat height: 480 mm
- Variants for artificial leather eligible for upholstery covering (see below)
- Load capacity: 200 kg
- Weight: 6 kg



### Model 404-01 chromium-plated

- Not movable
- 5-leg-stool, made of round pipe, Ø 520 mm, chromium-plated
- Seat plate made of laminated wood, Ø 350 mm
- Seat flat upholstered, thickness: 35 mm
- Seat height: 490 mm
- Variants for artificial leather eligible for upholstery covering (see below)
- Load capacity: 200 kg
- Weight: 4 kg

Design:

- 2 different variants for artificial leather eligible for upholstery covering:
  - Skai® Pandoria Plus with velvet-smooth surface
  - Skai® Toronto EN with innovative staynu-technology



# Useful devices for daily practice

## Swivel stools with height adjustment

### Swivel stools - Models 410-02, 411-02, 412-02, 414-02

- 5-arm polished aluminium base, Ø 560 mm
- **Seat:**
  - Seat height adjustable by gas spring; around going hand ring release under seat for height positioning
  - Seat base made of laminated wood
- **Load capacity:** 130 kg







411-02K, 411-02, 411-02L; in medium height



#### Design:

- 2 different variants for artificial leather eligible for upholstery covering:
  - Skai® Pandora Plus with velvet-smooth surface
  - Skai® Toronto EN with innovative staynu-technology



Mobility	Seat upholstery	Moulded upholstery, Ø 370 mm Thickness: 70 mm
<p><b>Movable</b> on 5 high-quality universal rollers, Ø 50 mm; for hard and soft floors</p> <ul style="list-style-type: none"> <li>• Cover made of zinc-die-casting, chromium-plated</li> <li>• Rollers according to legal regulations with brackets, only movable under load</li> </ul>	<p style="text-align: center;"><b>Model 410-02</b></p> <ul style="list-style-type: none"> <li>• Height adjustment: 490 mm - 680 mm</li> <li>• <b>Weight:</b> 6,5 kg</li> </ul> <div style="text-align: center;">  </div> <p><b>Variants:</b></p> <ul style="list-style-type: none"> <li>• <b>410-02K</b> Height adjustment: 390 mm - 520 mm <span style="font-size: 2em; vertical-align: middle;">K</span></li> <li>• <b>410-02L</b> Height adjustment: 540 mm - 780 mm <span style="font-size: 2em; vertical-align: middle;">L</span></li> </ul>	<p style="text-align: center;"><b>Model 411-02</b></p> <ul style="list-style-type: none"> <li>• Height adjustment: 530 mm - 720 mm</li> <li>• <b>Weight:</b> 8 kg</li> </ul> <div style="text-align: center;">  </div> <p><b>Variants:</b></p> <ul style="list-style-type: none"> <li>• <b>411-02K</b> Height adjustment: 420 mm - 550 mm <span style="font-size: 2em; vertical-align: middle;">K</span></li> <li>• <b>411-02L</b> Height adjustment: 580 mm - 820 mm <span style="font-size: 2em; vertical-align: middle;">L</span></li> </ul>
<p><b>Not movable</b> Fixed standing on floor glides with high bearing face (50 mm)</p>	<p style="text-align: center;"><b>Model 412-02</b></p> <ul style="list-style-type: none"> <li>• Height adjustment: 460 mm - 650 mm</li> <li>• <b>Weight:</b> 6,5 kg</li> </ul> <div style="text-align: center;">  </div>	<p style="text-align: center;"><b>Model 414-02</b></p> <ul style="list-style-type: none"> <li>• Height adjustment: 500 mm - 690 mm</li> <li>• <b>Weight:</b> 8 kg</li> </ul> <div style="text-align: center;">  </div>

## Swivel stools with height adjustment

### Swivel stool with foot operated height adjustment - Model 417-02

- 5-arm polished aluminium base
- **Mobility:**
  - Movable on 5 high-quality universal rollers, Ø 65 mm; for hard and soft floors
  - Cover made of zinc-die-casting, chromium-plated
  - Rollers according to legal regulations with brackets, only movable under load
- **Seat:**
  - Seat height adjustable by gas spring: 540 mm - 720 mm; Foot operated height adjustment by push plate at the lower frame
  - Seat base made of laminated wood
  - Moulded upholstery, Ø 370 mm, thickness 70 mm
- **Load capacity:** 130 kg
- **Weight:** 9 kg



### Swivel Stool - Model 420-02

- 5-arm polished aluminium base
- **Mobility:**
  - Movable on 5 high-quality universal rollers, Ø 65 mm; for hard and soft floors
  - Cover made of zinc-die-casting, chromium-plated
  - Rollers according to legal regulations with brackets, only movable under load
- With going around foot rest
- **Seat:**
  - Seat height adjustable by gas spring: 550 mm - 740 mm; with release hand-lever under seat
  - Seat base made of laminated wood
  - Moulded upholstery, Ø 420 mm, thickness 70 mm
  - Upholstered backrest, manual height adjustable
- **Load capacity:** 130 kg
- **Weight:** 12 kg



### Variants:



**420-02K**  
Height adjustment:  
440 mm - 570 mm

**K**



**420-02L**  
Height adjustment:  
600 mm - 840 mm

**L**

### Design:

- 2 different variants for artificial leather eligible for upholstery covering:
  - Skai® Pandoria Plus with velvet-smooth surface
  - Skai® Toronto EN with innovative staynu-technology







## Useful devices for daily practice

Steps made of square tube

### Technical data:

- Frame made of square tube, 35 x 35 x 3 mm
- Platforms: PUR-Foam with wooden basis and grained surface; colour: black
- Platform(s): 550 x 300 mm



Frame colours  
(left to right: -00, -03, -04)

**Model 319-00 white powder coated (RAL 9010)**

**Model 319-03 white aluminium powder coated (RAL 9006)**

**Model 319-04 grey aluminium powder coated (RAL 9007)**

### Measures:

- Platform: 550 x 300 mm
- Height: 150 mm
- Weight: 4,5 kg



Model 319-00

**Model 320-00 white powder coated (RAL 9010)**

**Model 320-03 white aluminium powder coated (RAL 9006)**

**Model 320-04 grey aluminium powder coated (RAL 9007)**

### Measures:

- Platform: 550 x 300 mm
- Height: 250 mm
- Weight: 5 kg



**Model 321-00 white powder coated (RAL 9010)**

**Model 321-03 white aluminium powder coated (RAL 9006)**

**Model 321-04 grey aluminium powder coated (RAL 9007)**

### Measures:

- Platforms: 550 x 300 mm
- Height Step 1: 230 mm
- Height Step 2: 460 mm
- Length overall: 560 mm
- Width overall: 300 mm
- Weight: 11 kg



Model 321-03

# Useful devices for daily practice

## Steps with handrail

**Model 323-00 white powder coated (RAL 9010)**

**Model 323-03 white aluminium powder coated (RAL 9006)**

**Model 323-04 grey aluminium powder coated (RAL 9007)**



### Equipment:

- With **one-sided handrail** for help and safety of your patients
- Without extra charge: Handrail right or left



### Technical data:

- Frame made of square tube 35 x 35 x 3 mm
- Platforms: PUR-Foam with wooden basis; colour: black
- Frame additionally with lateral support to prevent the risk of tipping over



### Measures:

- Height step 1: 230 mm | Height step 2: 460 mm | Platforms: 550 x 300 mm
- Height upper edge handrail: 1.350 mm | Width incl. lateral support: 800 mm | Depth: 600 mm



**Model 324-00 white powder coated (RAL 9010)**

**Model 324-03 white aluminium powder coated (RAL 9006)**

**Model 324-04 grey aluminium powder coated (RAL 9007)**



### Equipment:

- With **both-sided handrail** for help and safety of patients
- Without extra charge: Handrail right or left



### Technical data:

- Frame made of square tube 35 x 35 x 3 mm
- Platforms: PUR-Foam with wooden basis; colour: black
- Frame additionally with lateral support to prevent the risk of tipping over

### Measures:

- Height step 1: 230 mm
- Height step 2: 460 mm
- Height upper edge handrail: 1.350 mm
- Platforms: 550 x 300 mm
- Width incl. lateral support: 1.050 mm
- Depth: 600 mm

## Upholstery protection and care



### Cleaning kit for Skai products (Art. no. 050-1000)

for careful cleaning of upholstery

Consists of:


- Spray bottle with 200 ml special Skai cleansing agent
- 1 microfibre care cloth 40 x 40 cm



### Protection mat

Art. no. 1890

made of PVC, thickness: 0.8 mm,  
measures: 700 x 600 mm



K.H. Dewert GmbH  
Vollmestraße 7  
33649 Bielefeld  
Germany

Post office 14 40 28  
33634 Bielefeld  
Germany

Phone + 49 (0) 5 21 / 4 00 27-0  
Fax + 49 (0) 5 21 / 4 00 27-27

Internet: [www.khdewert.de](http://www.khdewert.de)  
E-mail: [info@khdewert.de](mailto:info@khdewert.de)  
YouTube: K.H. Dewert

04.2023\_02

USEFUL DEVICES FOR DAILY PRACTICE

