


Series 2000XL

3-sectioned versatile treatment tables with hydraulic height adjustment

… for examination, treatment, therapy …

Model 2008-00XL/H and 2009-00XL/H white powder coated (RAL 9010)
Model 2008-03XL/H and 2009-03XL/H white aluminium powder coated (RAL 9006)
Model 2008-04XL/H and 2009-04XL/H silver aluminium powder coated (RAL 9007)

▶Model 2808-..XL/H corresponds with model 2008-..XL/H, but with nose opening



German manufacture

All pictures:
Model 2008-00XL/H
with option:
Nose opening and
upholstery part to close
(1610 and 1620);
upholstery colour: 64-717 deep sea

Basic equipment:

- Table not movable, including height adjustment screws to compensate for uneven ground
- Foot part can be hinged vertical
- Pedals for foot-controlled height positioning on both sides of the table
- Dewert-hygiene-equipment



Upholstery:



Lying surface 3-sectioned:

- Head part stepless up- (+50°) and downwards (-45°) adjustable by gas spring
- Length head part: 450 mm
- Fixed middle part
- Length middle part: 1290 mm
- Foot part stepless up- (+30°) and downwards (-90°) adjustable by gas spring
- Length foot part: 450 mm
- Upholstery thickness: 60 mm

Technical data:

- Manual foot operated height adjustment: 540 mm to 960 mm
  by a hydraulic pump with a pressure capacity of 10000 N
- Width of lying surface: Model 2008-..XL/H: 700 mm | Model 2009-..XL/H: 800 mm
- Length of lying surface: min. 1790 mm - max. 2190 mm
- Load capacity: 250 kg; according to IEC 60601 tested with fourfold static load




Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology





Available as an option (see pricelist for additional options)

Mobility



Art. no. 1120N-10
Movable on comfort twin castors Ø 100 mm, each with total lock, basic height -15 mm



Art. no. 1180N
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm



Art. no. 1185
Movable on comfort twin castors, Ø 125 mm with central lock, basic height + 50 mm



Art. no. 1045N
Wheel-lifting/fixing mechanism with 4 foot levers for controlled lowering

Paper roll holder



Art. no. 1270K/F / 1280K/F
Paper roll holder, mounted at foot part (not protruding), for width 700 mm / 800 mm

Wall protection bumpers



Art. no. 1175 / 1176
4 wall protection bumpers Ø 150 mm, including mounting, for width 700 mm / 800 mm

Nose opening



Art. no. 1610
Nose opening, in head part

Art. no. 1620
Upholstery part to close

Upholstery options

Art. no. 1770 / 1780
Soft covering for pleasant resting comfort, for width 700 mm / 800 mm, thickness 80 mm

Art. no. 1772 / 1782
Upholstery medium-soft (cold foam cover), for width 700 mm / 800 mm, thickness 70 mm

Updated 06.2022_01
changes concerning the technical data are subject to change