

German manufacture

# Series 2000XL

## 2-sectioned versatile treatment tables with hydraulic height adjustment

**… for examination, treatment, therapy …**

**Model 2030-00XL/H and 2035-00XL/H white powder coated (RAL 9010)**
**Model 2030-03XL/H and 2035-03XL/H white aluminium powder coated (RAL 9006)**
**Model 2030-04XL/H and 2035-04XL/H silver aluminium powder coated (RAL 9007)**




All pictures:
Model 2035-00XL/H
with options: Movable on twin castors with central lock (1180N), paper roll holder (1280K);
upholstery colour:
64-717 deep sea

### Basic equipment:

- Table not movable, including height adjustment screws to compensate for uneven ground
- Long head part with half-round design
- Pedals for foot-controlled height positioning on both sides of the table
- Dewert-hygiene-equipment




### Upholstery:

Lying surface 2-sectioned:

- Head part stepless adjustable from -20° to +70° by gas spring (sitting position)
- Length head part: 750 mm
- Fixed foot part
- Length foot part: 1.200 mm
- Upholstery thickness: 60 mm

### Technical data:


- Manual foot operated height adjustment: 540 mm to 960 mm
  by a hydraulic pump with a pressure capacity of 10000 N
- Width of lying surface: Model 2030-..XL/H: 700 mm | Model 2035-..XL/H: 800 mm
- Length of lying surface: 1950 mm
- Load capacity: 250 kg; according to IEC 60601 tested with fourfold static load

kg

### Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



vertical
adjustment

## Available as an option (see pricelist for additional options)

### Mobility



Art. no. 1120N-10
Movable on comfort twin castors Ø 100 mm, each with total lock, basic height -15 mm



Art. no. 1180N
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm



Art. no. 1185
Movable on comfort twin castors, Ø 125 mm with central lock, basic height + 50 mm



Art. no. 1045N
Wheel-lifting/fixing mechanism with 4 foot levers for controlled lowering

### Paper roll holder



Art. no. 1270K / 1280K
Paper roll holder, mounted under head part (not protruding), for width 700 mm / 800 mm



Art. no. 1270F / 1280F
Paper roll holder, mounted at foot part (protruding: 70 mm), for width 700 mm / 800 mm

### Wall protection bumpers



Art. no. 1175 / 1176
4 wall protection bumpers Ø 150 mm, including mounting, for width 700 mm / 800 mm

### Side guards



Art. no. 1080-60
1 side guard lateral lowerable, chromium-plated, completely made of steel tube, length: 820 mm

### Nose opening



Art. no. 1610
Nose opening, in head part

Art. no. 1620
Upholstery part to close

### Upholstery options

Art. no. 1770 / 1780
Soft covering for pleasant resting comfort, for width 700 mm / 800 mm, thickness 80 mm

Art. no. 1772 / 1782
Upholstery medium-soft (cold foam cover), for width 700 mm / 800 mm, thickness 70 mm

Updated 06.2022_01
changes concerning the technical data are subject to change