


Series 2000XL

Short tables with hydraulic height adjustment

… for examination, treatment, therapy …

Model 2100-00XL/H, 2105-00XL/H, 2110-00XL/H white powder coated (RAL 9010)
Model 2100-03XL/H, 2105-03XL/H, 2110-03XL/H white aluminium powder coated (RAL 9006)
Model 2100-04XL/H, 2105-04XL/H, 2110-04XL/H silver aluminium powder coated (RAL 9007)



German manufacture



Both pictures:
Model 2100-00XL/H
with option: Movable on twin castors with central lock (1180N);
upholstery colour: 64-717 deep sea

Basic equipment:

- Table not movable, including height adjustment screws to compensate for uneven ground
- Pedals for foot-controlled height positioning on both sides of the table
- Dewert-hygiene-equipment

Upholstery:

Lying surface 1-sectioned:

- Length: 1400 mm (other measures on request)
- Upholstery thickness: 60 mm



Table in highest position

Technical data:

- Manual foot operated height adjustment: 50 mm to 920 mm
  by a hydraulic pump with a pressure capacity of 10000 N
- Width of lying surface: Model 2100-..XL/H: 650 mm | Model 2105-..XL/H: 800 mm | Model 2110-..XL/H: 1.000 mm
- GesamtLength of lying surface: 1400 mm
- Load capacity: 250 kg; according to IEC 60601 tested with fourfold static load




Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology





Available as an option (see pricelist for additional options)

Mobility



Art. no. 1120N-10
Movable on comfort twin castors Ø 100 mm, each with total lock, basic height -15 mm



Art. no. 1180N
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm



Art. no. 1185
Movable on comfort twin castors, Ø 125 mm with central lock, basic height + 50 mm



Art. no. 1045N
Wheel-lifting/fixing mechanism with 4 foot levers for controlled lowering

DEWERT
Präzisions-Stahlrohrmöbel

Updated 06.2022_01
changes concerning the technical data are subject to change