Series DFH

2-sectioned table with fixed height made of round tube
Length head part: 550 mm

Model 106-01 and 107-01 chromium-plated
Model 106-00 and 107-00 white powder coated (RAL 9010)
Model 106-03 and 107-03 white aluminium powder coated (RAL 9006)
Model 106-04 and 107-04 grey aluminium powder coated (RAL 9007)


German manufacture



Both pictures:
Model 106-04 with option:
Paper roll holder (1210);
upholstery colour: 64-717 deep sea

Basic equipment:

- Frame made of round tube
- Table not movable, including height adjustment screws to compensate for uneven ground
- Head part positive adjustable
- Dewert-hygiene-equipment



Upholstery:

Lying surface 2-sectioned:

- Head part upwards (+30°) adjustable by 2 metal rachets
- Length head part: 550 mm
- Fixed foot part
- Length foot part: 1400 mm
- Upholstery thickness: 80 mm



Technical data:

- Length overall: 1950 mm
- Width overall: Model 106-..: 650 mm | Model 107-..: 800 mm
- Height: 650 mm
- Load capacity: 250 kg kg

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology

Available as an option (see pricelist for additional options)

Mobility



Art. no. 1120-10
Movable on 4 castors, Ø 100 mm, each with total lock



Art. no. 1140-11
Movable on 4 castors, Ø 125 mm, each with total lock

Longitudinal brace



Art. no. 1310
Longitudinal brace for additional stability, necessary in case of mobility

Paper roll holder



Art. no. 1210 / 1220
Paper roll holder, for width 650 mm / 800 mm

Upholstery options

Art. no. 1765 / 1780
Soft covering for pleasant resting comfort, for width 650 mm / 800 mm, thickness 80 mm

Art. no. 1762 / 1782
Upholstery medium-soft (cold foam cover), for width 650 mm / 800 mm, thickness 70 mm

Head part options



Art. no. 1655
Head part stepless adjustable by gas spring



Art. no. 1662
Additional second head part of the same size, adjustable by 2 metal rachets



Art. no. 1610
Nose opening

Art. no. 1620
Upholstery part to close

Updated 06.2022_01
changes concerning the technical data are subject to change