

Series DFH

Short tables for examination and treatment
3-sectioned table with fixed height (square tube)

Model 113-00 and 114-00 white powder coated (RAL 9010)
Model 113-01 and 114-01 chromium-plated
Model 113-03 and 114-03 white aluminium powder coated (RAL 9006)
Model 113-04 and 114-04 grey aluminium powder coated (RAL 9007)

German manufature



All pictures:
Model 113-04 with basic equipment;
upholstery colour: 64-714 olive green



Basic equipment:

- Foot part can be hinged vertical
- Table not movable, including height adjustment screws to compensate for uneven ground
- Dewert-hygiene-equipment

Upholstery:

Lying surface 3-sectioned:

- Head part upwards (+30°) adjustable by 2 metal rachets
- Length head part: 550 mm
- Fixed middle part
- Length middle part: 920 mm
- Foot part stepless adjustable (-90° - +10°) by gas spring (can be hinged vertical)
- Length foot part: 480 mm
- Upholstery thickness: 80 mm



Technical data:

- Length overall: 1950 mm
- Width overall: Model 113-..: 650 mm | Model 114-..: 800 mm
- Height: 800 mm
- Load capacity: 250 kg kg

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology

Available as an option (see pricelist for additional options)

DEWERT Präzisions-Stahlrohrmöbel

Mobility



Art. no. 1120-10
Movable on 4 castors, Ø 100 mm, each with total lock

Longitudinal brace



Art. no. 1310
Longitudinal brace for additional stability, necessary in case of mobility

Paper roll holder



Art. no. 1210 / 1220
Paper roll holder, for width 650 mm / 800 mm

Upholstery options

Art. no. 1765 / 1780
Soft covering for pleasant resting comfort, for width 650 mm / 800 mm, thickness 80 mm

Art. no. 1762 / 1782
Upholstery medium-soft (cold foam cover), for width 650 mm / 800 mm, thickness 70 mm

Head part options



Art. no. 1655-20
Head part stepless adjustable by gas spring (0° to +55°)



Art. no. 1610
Nose opening

Art. no. 1620
Upholstery part to close

Side guards



Art. no. 1088
2 side guards lateral lowerable, RAL 9006 powder coated, length: 820 mm



Art. no. 1040
1 side guard tiltable, chromium-plated, completely made of steel tube, length: 700 mm



Art. no. 1005
1 side guard to fit on side rail, with attachment clamp, chromium-plated, length: 800 mm,

Add-on elements



1 piece of side rail, length: 350 mm

Art. no. 1420
1 piece of side rail, length: 1300 mm



Art. no. 1430
1 attachment clamp for side rail, suitable for pipe: Ø 18 mm



Art. no. 1440
1 arm rest for side rail, without attachment clamp



Art. no. 1450
Infusion pole, without attachment clamp

Updated 06.2022_01
changes concerning the technical data are subject to change