

Series DFH

3-sectioned table with fixed height
Length head part: 550 mm

Model 120-00 white powder coated (RAL 9010)
Model 120-01 chromium-plated
Model 120-03 white aluminium powder coated (RAL 9006)
Model 120-04 grey aluminium powder coated (RAL 9007)

German manufacture



Both pictures:
Model 120-00
with option:
Paper roll holder (1210);
upholstery colour: 64-701 stone blue


Basic equipment:

- Frame made of round tube
- Table not movable, including height adjustment screws to compensate for uneven ground
- Longitudinal brace for additional stability
- Head part positive adjustable
- Foot part positive adjustable
- Dewert-hygiene-equipment

Upholstery:

Lying surface 3-sectioned:

- Head part upwards (+30°) adjustable by 2 metal rachets
- Length head part: 550 mm
- Fixed middle part
- Length middle part: 450 mm
- Foot part stepless adjustable by gas spring (0° to +28°)
- Length foot part: 950 mm
- Upholstery thickness: 60 mm



Technical data:

- Length overall: 1950 mm
- Width overall: 650 mm
- Height: 650 mm
- Load capacity: 250 kg

kg

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



Available as an option (see pricelist for additional options)

Mobility



Art. no. 1120-10
Movable on 4 castors, Ø 100 mm, each with total lock



Art. no. 1140-10
Movable on 4 castors, Ø 125 mm, each with total lock

Paper roll holder



Art. no. 1210
Paper roll holder, for width 650 mm / 800 mm

Upholstery options

Art. no. 1765
Soft covering for pleasant resting comfort, for width 650 mm, thickness 80 mm

Art. no. 1762
Upholstery medium-soft (cold foam cover), for width 650 mm, thickness 70 mm

Head part options



Art. no. 1655
Head part stepless adjustable by gas spring



Art. no. 1610
Nose opening

Art. no. 1620
Upholstery part to close

Updated 06.2022_01
changes concerning the technical data are subject to change