# Series DFH

## Short tables for examination and treatment
### 3-sectioned table with fixed height (round tube)

**Model 121-01 chromium-plated**
**Model 121-00 white powder coated (RAL 9010)**
**Model 121-03 white aluminium powder coated (RAL 9006)**
**Model 121-04 grey aluminium powder coated (RAL 9007)**

German manufacture



All pictures:
Model 121-01
with option:
Paper roll holder (1210);
upholstery colour: 64-714 olive green




### Basic equipment:

- Frame made of round tube
- Longitudinal brace for additional stability
- Foot part can be hinged vertical
- Table not movable, including height adjustment screws to compensate for uneven ground
- Dewert-hygiene-equipment

### Upholstery:




Lying surface 3-sectioned:

- Head part upwards (+30°) adjustable by 2 metal rachets
- Length head part: 550 mm
- Fixed middle part
- Length middle part: 900 mm
- Foot part stepless adjustabe (-90° to +20°) by gas spring (can be hinged vertical)
- Length foot part: 500 mm
- Upholstery thickness: 80 mm

### Technical data:

- Length overall: min. 1550 mm, max. 1970 mm
- Width overall: 650 mm
- Height: 800 mm
- Load capacity: 250 kg (foot part vertical downwards)

kg

### Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology

## Available as an option (see pricelist for additional options)

### Mobility



Art. no. 1120-10
Movable on 4 castors, Ø 100 mm, each with total lock

### Head part options



Art. no. 1655
Head part stepless adjustable by gas spring



Art. no. 1610
Nose opening

Art. no. 1620
Upholstery part to close

### Paper roll holder



Art. no. 1210
Paper roll holder,
for width 650 mm / 800 mm

### Upholstery options

Art. no. 1765
Soft covering for pleasant resting comfort, for width 650 mm, thickness 80 mm

Art. no. 1762
Upholstery medium-soft (cold foam cover), for width 650 mm, thickness 70 mm

Updated 06.2022_01
changes concerning the technical data are subject to change