Series DFH

Short tables for examination and treatment

3-sectioned table with fixed height (round tube)

Model 122-00 white powder coated (RAL 9010)
Model 122-01 chromium-plated
Model 122-03 white aluminium powder coated (RAL 9006)
Model 122-04 grey aluminium powder coated (RAL 9007)

German manufacture





All pictures.:
Model 122-00
with options:
Movable on castors, each with lock (1120-10),
Paper roll holder (1265K);
upholstery colour: 64-708 cherry

Basic equipment:

- Frame made of round tube
- Head part and foot part can be hinged vertical
- Table not movable, including height adjustment screws to compensate for uneven ground
- Longitudinal brace for additional stability
- Dewert-hygiene-equipment



Upholstery:

Lying surface 3-sectioned:

- Head part stepless adjustabe (-90° to +25°) by gas spring (can be hinged vertical)
- Length head part: 450 mm
- Fixed middle part
- Length middle part: 1120 mm
- Foot part stepless adjustabe (-90° to +25°) by gas spring (can be hinged vertical)
- Length foot part: 450 mm
- Upholstery thickness: 80 mm

Technical data:

- Length overall: min. 1300 mm, max. 2020 mm
- Width overall: 650 mm
- Height: 800 mm
- Load capacity: 250 kg (head part and foot part vertical downwards)

kg

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology

Available as an option (see pricelist for additional options)

Mobility



Art. no. 1120-10
Movable on 4 castors,
Ø 100 mm, each with total lock


Head part options



Art. no. 1610
Nose opening

Art. no. 1620
Upholstery part to close

Paper roll holder



Art. no. 1265K
Paper roll holder, for width 650 mm

Upholstery options

Art. no. 1765
Soft covering for pleasant resting comfort, for width 650 mm, thickness 80 mm

Art. no. 1762
Upholstery medium-soft (cold foam cover), for width 650 mm, thickness 70 mm

DEWERT Präzisions-Stahlrohrmöbel

Updated 06.2022_01
changes concerning the technical data are subject to change