Series DFH

2-sectioned table with fixed height made of square tube
Length head part: 550 mm

Model 27-00 and 37-00 white powder coated (RAL 9010)
Model 27-03 and 37-03 white aluminium powder coated (RAL 9006)
Model 27-04 and 37-04 grey aluminium powder coated (RAL 9007)

German manufacture



Both pictures:
Model 27-04 with option:
Paper roll holder (1210);
upholstery colour: 64-715 anthracite

Basic equipment:

- Head part positive and negative adjustable
- Longitudinal brace for additional stability
- Frame welded as one unit
- Table not movable, including height adjustment screws to compensate for uneven ground
- Dewert-hygiene-equipment

Upholstery:

Lying surface 2-sectioned:

- Head part stepless adjustable by gas spring (-25° to +30°)
- Length head part: 550 mm
- Fixed foot part
- Length foot part: 1400 mm
- Upholstery thickness: 60 mm



Technical data:

- Length overall: 1950 mm
- Width overall: Model 27-..: 650 mm | Model 37-..: 800 mm
- Height: 650 mm
- Load capacity: 250 kg

kg

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



Available as an option (see pricelist for additional options)

Mobility



Art. no. 1120-10
Movable on 4 castors, Ø 100 mm, each with total lock



Art. no. 1140-10
Movable on 4 castors, Ø 125 mm, each with total lock



Art. no. 1180N
Movable on twin castors, Ø 100 mm with central lock



Art. no. 1185
Movable on twin castors, Ø 125 mm with central lock

DEWERT Präzisions-Stahlrohrmöbel

Head part options



Art. no. 1610
Nose opening

Art. no. 1620
Upholstery part to close

Paper roll holder



Art. no. 1210 / 1220
Paper roll holder, for width 650 mm / 800 mm

Upholstery options

Art. no. 1765 / 1780
Soft covering for pleasant resting comfort, for width 650 mm / 800 mm, thickness 80 mm

Art. no. 1762 / 1782
Upholstery medium-soft (cold foam cover), for width 650 mm / 800 mm, thickness 70 mm

Side guards



Art. no. 1088
2 side guards lateral lowerable, RAL 9006 powder coated, length: 820 mm



Art. no. 1040
1 side guard tiltable, chromium-plated, completely made of steel tube, length: 700 mm



Art. no. 1005
1 side guard to fit on side rail, with attachment clamp, chromium-plated, length: 800 mm,

Add-on elements



Art. no. 1410
1 piece of side rail, length: 350 mm

Art. no. 1420
1 piece of side rail, length: 1300 mm



Art. no. 1430
1 attachment clamp for side rail, suitable for pipe: Ø 18 mm



Art. no. 1440
1 arm rest for side rail, without attachment clamp



Art. no. 1450
Infusion pole, without attachment clamp

Updated 06.2022_01
changes concerning the technical data are subject to change