

Series DFH

3-sectioned table with fixed height made of square tube

Length head part: 750 mm

Model 30-00 white powder coated (RAL 9010)
Model 30-03 white aluminium powder coated (RAL 9006)
Model 30-04 grey aluminium powder coated (RAL 9007)

NEW

German manufacture



Both pictures:
Model 30-00
basic equipment;
upholstery colour: 64-717 deep sea

Basic equipment:

- Table not movable, including height adjustment screws to compensate for uneven ground
- Longitudinal brace for additional stability
- Head part and foot part positive and negative adjustable
- Possibility to pin on a paper roll holder on head side and foot side
- Dewert-hygiene-equipment

Upholstery:

Lying surface 3-sectioned:

- Head part stepless adjustable (0° to +80°) by gas spring (sitting position)
- Length head part: 750 mm
- Fixed middle part
- Length middle part: 300 mm
- Foot part stepless adjustable by gas spring (0° to +30°)
- Length foot part: 900 mm
- Upholstery thickness: 60 mm



Technical data:

- Length overall: 1950 mm
- Width overall: 650 mm, optional 800 mm
- Height: 680 mm
- Load capacity: 250 kg

kg

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



Available as an option (see pricelist for additional options)

DEWERT Präzisions-Stahlrohrmöbel

Mobility



Art. no. 1120-10
Movable on 4 castors, Ø 100 mm, each with total lock



Art. no. 1140-10
Movable on 4 castors, Ø 125 mm, each with total lock



Art. no. 1180N
Movable on twin castors, Ø 100 mm with central lock



Art. no. 1185
Movable on twin castors, Ø 125 mm with central lock


Head part options



Art. no. 1610
Nose opening

Art. no. 1620
Upholstery part to close

Side guards



Art. no. 1088
2 side guards lateral lowerable, RAL 9006 powder coated, length: 820 mm



Art. no. 1040
1 side guard tiltable, chromium-plated, completely made of steel tube, length: 700 mm



Art. no. 1005
1 side guard to fit on side rail, with attachment clamp, chromium-plated, length: 800 mm,

Paper roll holder



Art. no. 1210 / 1220
Paper roll holder,
for width 650 mm / 800 mm

Upholstery options

Art. no. 30-800:
Width of upholstery 800 mm

Art. no. 1765 / 1780
Soft covering for pleasant resting comfort, for width 650 mm / 800 mm, thickness 80 mm

Art. no. 1762 / 1782
Upholstery medium-soft (cold foam cover), for width 650 mm / 800 mm, thickness 70 mm



Art. no. 1850
Storage basket with mounting, measures:
L. 530 x W. 370 x H. 80 mm

Add-on elements



Art. no. 1410
1 piece of side rail, length: 350 mm

Art. no. 1420
1 piece of side rail, length: 1300 mm



Art. no. 1430
1 attachment clamp for side rail, suitable for pipe: Ø 18 mm



Art. no. 1440
1 arm rest for side rail, without attachment clamp



Art. no. 1450
Infusion pole, without attachment clamp

Updated 06.2022_01
changes concerning the technical data are subject to change