


# Series DX1

## 3-sectioned versatile treatment tables with electric motorized height adjustment

**...for examination, treatment, therapy…**

**Model 1050-00E white powder coated (RAL 9010)**
**Model 1050-03E white aluminium powder coated (RAL 9006)**
**Model 1050-04E grey aluminium powder coated (RAL 9007)**




German manufacture



Both pictures:
Model 1050-00E with options:
Movable on twin castors with central lock (046), loss of nose opening (1615-EF);
upholstery colour: 64-717 deep sea

### Basic equipment:

- Table not movable, including height adjustment screws to compensate for uneven ground
- Electromotive drive control by foot switch (option: hand switch)
- Integrated safety function for automatic blocking of electromotive drive
- Freewheeling downwards, automatic declutch in case of resistance
- Dewert-hygiene-equipment

### Upholstery:

Lying surface 3-sectioned:

- Head part stepless up- (+40°) and downwards (-30°) adjustable by gas spring
- Length head part: 550 mm, integrated nose opening and upholstery part to close
- Fixed middle part
- Length middle part: 480 mm
- Foot part stepless adjustable by gas spring (0° to +50°)
- Length foot part: 920 mm
- Upholstery thickness: 65 mm




### Technical data:

- Electronical height adjustment: 500 - 940 mm in 22 sec. by AC-motor with a pressure capacity of 10000 N
- Width overall: 700 mm
- Length overall: 1950 mm
- Load capacity: 225 kg, according to IEC 60601 tested with fourfold static load
  Load capacity head part: 90 kg


kg

### Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology






## Available as an option (see pricelist for additional options)

DEWERT Präzisions-Stahlrohrmöbel

### Mobility



Art. no. 040N
Movable on comfort twin castors Ø 100 mm, each with total lock, basic height -15 mm



Art. no. 046
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm



Art. no. B1185
Movable on comfort twin castors, Ø 125 mm with central lock, basic height + 50 mm



Art. no. 045N
Wheel-lifting/fixing mechanism with 4 foot levers for controlled lowering

### Additional motor control options



Art. no. 090E
Hand switch (low voltage), with magnet on the back to fasten



Art. no. 091E
Foot switch (low voltage)



Art. no. 092E
1 foot switch fixation mounted at underframe



Art. no. 094E
Foot switch rail, lateral, on both sides, security adjustment:
Press the switch rail down = the table is lifted
Lift the switch rail up = the table is lowered



Art. no. 096E
Foot switch rail, surround, security adjustment:
Press the switch rail down = the table is lifted
Lift the switch rail up = the table is lowered

### Additional safety options



Art. no. B1997E
Foot switch cover, made of brushed VA-steel

### Paper roll holder



Art. no. 010K/015K
Paper roll holder, mounted under head part (not protruding), for width 700 mm / 800 mm




Art. no. 010K/F / 015K/F
Paper roll holder, mounted at foot part (not protruding), for width 700 mm / 800 mm

### Wall protection bumpers



Art. no. B1175 / B1176
4 wall protection bumpers Ø 150 mm, including mounting, for width 700 mm / 800 mm

### Head part options



Art. no. 1615-EF
Loss of nose opening and upholstery part to close



Art. no. 056
Head part 3-sectioned, arm rests downwards adjustable and removable
(Art. no. 010K/015K not possible)

### Side guards



Art. no. 1008
2 side guards lateral lowerable, RAL 9006 powder coated, made of steel tube, length: 820 mm

### Upholstery options

Art. no. 007: Width 800 mm

Art. no. B1770 / B1788:
Soft covering for pleasant resting comfort, for width 700 mm/ 800 mm, thickness 80 mm

Art. no. B1772 / B1782:
Upholstery medium-soft (cold foam cover), for width 700 mm / 800 mm, thickness 70 mm

Updated 2022.1_01
Technical data are subject to change