


Serie DX1

3-sectioned versatile treatment tables with hydraulic height adjustment

...for examination, treatment, therapy…

Model 1050-00/H white powder coated (RAL 9010)
Model 1050-03/H white aluminium powder coated (RAL 9006)
Model 1050-04/H grey aluminium powder coated (RAL 9007)



German manufacture



Both pictures:
Model 1050-03/H with options:
Movable on comfort twin castors, each with lock (040N), without nose opening (1615-EF); upholstery colour: 64-7716 black

Basic equipment:

- Table not movable, including height adjustment screws to compensate for uneven ground
- Pedals for foot-controlled height positioning on both sides of the table
- Dewert-hygiene-equipment

Upholstery:

Lying surface 3-sectioned:

- Head part stepless up- (+40°) and downwards (-30°) adjustable by gas spring
- Length head part: 550 mm, integrated nose opening and upholstery part to close
- Fixed middle part
- Length middle part: 480 mm
- Foot part stepless adjustable by gas spring (0° to +50°)
- Length foot part: 920 mm
- Upholstery thickness: 65 mm



Technical data:

- Manual foot operated height adjustment: 490 - 940 mm
  by a hydraulic pump with a pressure capacity of 10000 N
- Width overall: 700 mm
- Length overall: 1950 mm
- Load capacity: 225 kg, according to IEC 60601 tested with fourfold static load
  Load capacity head part: 90 kg


kg

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology





Available as an option (see pricelist for additional options)

DEWERT Präzisions-Stahlrohrmöbel

Mobility



Art. no. 040N
Movable on comfort twin castors Ø 100 mm, each with total lock, basic height -15 mm



Art. no. 046
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm



Art. no. B1185
Movable on comfort twin castors, Ø 125 mm with central lock, basic height + 50 mm



Art. no. 045N
Wheel-lifting/fixing mechanism with 4 foot levers for controlled lowering

Paper roll holder



Art. no. 010K/015K
Paper roll holder, mounted under head part (not protruding), for width 700 mm / 800 mm



Art. no. 010K/F / 015K/F
Paper roll holder, mounted at foot part (not protruding), for width 700 mm / 800 mm

Side guards



Art. no. 1008
2 side guards lateral lowerable, RAL 9006 powder coated, made of steel tube, length: 820 mm

Wall protection bumpers



Art. no. B1175 / B1176
4 wall protection bumpers Ø 150 mm, including mounting, for width 700 mm / 800 mm

Head part options



Art. no. 056
Head part 3-sectioned, arm rests downwards adjustable from -25° to +45° and removable (Art. no. 010K/015K not possible)
measures:
Head part: L. 550 mm x W. 380 mm
Arm rests: L. 500 mm x W. 160 mm



Art. no. 1615-EF
Loss of nose opening and upholstery part to close

Upholstery options

Art. no. 007
Width 800 mm

Art. no. B1770 / B1788
Soft covering for pleasant resting comfort, for width 700 mm/ 800 mm, thickness 80 mm

Art. no. B1772 / B1782
Upholstery medium-soft (cold foam cover), for width 700 mm / 800 mm, thickness 70 mm



Updated 2022.1_01
Technical data are subject to change