


German manufacture

Series 2800XL

Gymnastic table for therapy acc. to Bobath and Vojta with electric motorized height adjustment

Model 2800-00XL/H* white powder coated (RAL 9010)
Model 2800-03XL/H* white aluminium powder coated (RAL 9006)
Model 2800-04XL/H* grey aluminium powder coated (RAL 9007)

* Exact name of model depending on width of lying surface (see technical data)



without picture



Model 2800-00XLE with option:
Movable on twin castors with central lock (046);
upholstery colour: 64-717 deep sea

Basic equipment:

- Lying surface 1-sectioned
- Table not movable, including height adjustment screws to compensate for uneven ground
- Pedals for foot-controlled height positioning on both sides of the table
- Dewert-hygiene-equipment



Upholstery:

Lying surface 1-sectioned:

- Length: 1950 mm
- Optional: Head part
- Thickness of upholstery: 60 mm



Model 2804-00XLE with options:
Wheel-lifting/fixing mechanism (045N), surrounding foot control (095E), head part stepless adjustable by gas spring (063);
upholstery colour: 64-707 chestnut

Technical data:

- Manual foot operated height adjustment: 500 mm to 920 mm by a hydraulic pump, pressure capacity of 10000 N
- Width of lying surface:
  Model 2800XL/H: 1000 mm | **Model 2802XL/H:** 1100 mm | **Model 2804XL/H:** 1200 mm
- Length of lying surface: 1950 mm
- Load capacity: 250 kg; according to IEC 60601 tested with fourfold static load


kg

Design:

2 different variants for artificial leather eligible for upholstery covering:


- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology





Available as an option (see pricelist for additional options)

Mobility



Art. no. 040N
Movable on comfort twin castors Ø 100 mm, each with total lock, basic height -15 mm



Art. no. 046
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm



Art. no. 045N
Wheel-lifting/fixing mechanism with 4 foot levers for controlled lowering

Paper roll holder



Art. no. 016K
Paper roll holder, mounted under head part (not protruding), for width 1000 mm

Head part options



Art. no. 063
Head part, upwards and downwards adjustable by gas spring