


German manufacture

Series 2800XL

3-sectioned therapy table with hydraulic height adjustment

Model 2820-00XL/H white powder coated (RAL 9010)
Model 2820-03XL/H white aluminium powder coated (RAL 9006)
Model 2820-04XL/H grey aluminium powder coated (RAL 9007)









All pictures:
Model 2820-04XL/H
with options:
Movable on twin castors with central lock (046), paper roll holder (010K/F);
upholstery colour: 64-713 corn

Basic equipment:

- Lying surface 3-sectioned
- Head part can be hinged vertical
- Foot part stepless upwards adjustable to sitting position (+80°)
- Table not movable, including height adjustment screws to compensate for uneven ground
- Pedals for foot-controlled height positioning on both sides of the table
- Dewert-hygiene-equipment

Upholstery:

Lying surface 3-sectioned:

- Head part stepless up- (+30°) and downwards (-90°) adjustable by gas spring
- Length head part: 450 mm, with nose opening and upholstery part to close
- Fixed middle part
- Length middle part: 530 mm
- Foot part stepless upwards (+80°) adjustable to sitting position by gas spring
- Length foot part: 1020 mm
- Thickness of upholstery: 60 mm



Technical data:

- Manual foot operated height adjustment: 530 mm to 950 mm by a hydraulic pump, pressure capacity of 10000 N
- Width of lying surface: 700 mm
- Length of lying surface: 2000 mm
- Load capacity: 250 kg; according to IEC 60601 tested with fourfold static load


kg

Design:

2 different variants for artificial leather eligible for upholstery covering:


- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology





Available as an option (see pricelist for additional options)

Mobility



Art. no. 040N
Movable on comfort twin castors Ø 100 mm, each with total lock, basic height -15 mm



Art. no. 046
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm



Art. no. B1185
Movable on comfort twin castors, Ø 125 mm with central lock, basic height + 50 mm



Art. no. 045N
Wheel-lifting/fixing mechanism with 4 foot levers for controlled lowering

Paper roll holder



Art. no. 010K
Paper roll holder, mounted under head part (not protruding), for width 700 mm



Art. no. 010K/F
Paper roll holder, mounted at foot part (not protruding), for width 700 mm

Upholstery options

Art. no. 007
Width of upholstery 800 mm

Art. no. B1770
Soft covering for pleasant resting comfort, for width 700 mm, thickness 80 mm

Art. no. B1772
Upholstery medium-soft (cold foam cover), for width 700 mm, thickness 70 mm

Add-on elements



Art. no. B1410
1 piece of side rail, length: 350 mm




Art. no. B1430
1 attachment clamp for side rail,
suitable for pipe: Ø 18 mm

Arm rests



Art. no. B1440
1 arm rest for side rail, without attachment clamp, swivelling and removable, additional revolving and tiltable

Nose opening



Art. no. 1615-EF
Loss of nose opening and upholstery part to close