


German manufacture

Series 2800XL

6-sectioned therapy table with hydraulic height adjustment

Model 2880-00XL/H white powder coated (RAL 9010)
Model 2880-03XL/H white aluminium powder coated (RAL 9006)
Model 2880-04XL/H grey aluminium powder coated (RAL 9007)



without picture



Model 2880-00E (eletric motorized height adjustment) with options: Wheel-lifting/fixing mechanism (045N), surrounding foot control (095E); upholstery colour: 64-702 white

Basic equipment:

- Lying surface 6-sectioned
- Middle part adjustable to roof position, Trendelenburg position possible
- Foot part 4-sectioned, split transversal and longitudinal
- Table not movable, including height adjustment screws to compensate for uneven ground
- Pedals for foot-controlled height positioning on both sides of the table
- Dewert-hygiene-equipment

Upholstery:

Lying surface 6-sectioned:

- Head part stepless up- (+45°) and downwards (-40°) adjustable by gas spring
- Length head part: 540 mm; with nose opening and upholstery part to close
- Middle part upwards manually adjustable to roof position, supported by 2 gas springs, due to 2 guide rods which are fixed by screwing a clamping lever, this adjustment system is without clearance
- Length middle part: 530 mm
- Foot part 4-sectioned, split transversal and longitudinal
- Every foot section stepless adjustable (+60°) by a self-fixing system of 2 guide rods, without clearance (foot sections independently adjustable)
- Lower foot sections additionally stepless adjustable (-65°) by gas spring
- Length foot part: Upper foot section: 405 mm | Lower foot section: 485 mm
- Thickness of upholstery: 60 mm

Technical data:

- Manual foot operated height adjustment: 590 mm to 1010 mm by a hydraulic pump, pressure capacity of 10000 N
- Width of lying surface: 700 mm
- Length of lying surface: 1950 mm
- Load capacity: 200 kg; according to IEC 60601 tested with fourfold static load


kg

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology





Available as an option (see pricelist for additional options)

Mobility



Art. no. 040N
Movable on comfort twin castors Ø 100 mm, each with total lock, basic height -15 mm



Art. no. 046
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm



Art. no. 045N
Wheel-lifting/fixing mechanism with 4 foot levers for controlled lowering

DEWERT Präzisions-Stahlrohrmöbel

Paper roll holder



Art. no. 010K
Paper roll holder, mounted under head part (not protruding), for width 700 mm, not in case of 3-sectioned head part (Art. no. 056)

Head part options



Art. no. 056
Head part 3-sectioned, arm rests downwards adjustable from -25° to +45° and removable
(art. no. 010K not possible)
measures:
Head part: L. 540 mm x W. 370 mm
Arm rests: L. 500 mm x W. 160 mm

Nose opening



Art. no. 1615-EF
Loss of nose opening and upholstery part to close

Arm rests

without picture

Art. no. B088
1 arm rest, fixed, swivelling and removable, additional revolving and tiltable

Upholstery options

Art. no. 009
Width of upholstery 800 mm

Art. no. B1772
Upholstery medium-soft (cold foam cover), for width 700 mm, thickness 70 mm

Updated 04.2023_01
changes concerning technical data are subject to change