


German manufacture

Series 2800XL

3-sectioned therapy table with hydraulic height adjustment

Model 2890-00XL/H white powder coated (RAL 9010)
Model 2890-03XL/H white aluminium powder coated (RAL 9006)
Model 2890-04XL/H grey aluminium powder coated (RAL 9007)











All pictures:
Model 2890-00XL/H
with options:
Movable on twin castors with central lock (046), paper roll holder (010F/K);
upholstery colour:
64-704 birch

Basic equipment:

- Lying surface 3-sectioned
- Multifunctional and independent adjustment of upholstery sections
- Possibilities of positioning the surface parts: Roof shape, Trendelenburg, reverse Trendelenburg
- Table not movable, including height adjustment screws to compensate for uneven ground
- Pedals for foot-controlled height positioning on both sides of the table
- Dewert-hygiene-equipment



Upholstery:

Lying surface 3-sectioned:

- Head part stepless up- (+40°) and downwards (-35°) adjustable by gas spring
- Length head part: 420 mm; with nose opening and upholstery part to close
- Middle part stepless up- (+45°) and downwards (-15°) adjustable by gas spring
- Length middle part: 550 mm
- Foot part stepless up- (+50°) and downwards (-15°) adjustable by gas spring
- Length foot part: 980 mm
- Thickness of upholstery: 60 mm

Technical data:

- Manual foot operated height adjustment: 600 mm to 1020 mm by a hydraulic pump, pressure capacity of 10000 N
- Width of lying surface: 700 mm
- Length of lying surface: 1950 mm
- Load capacity: 200 kg; according to IEC 60601 tested with fourfold static load


kg

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology





Available as an option (see pricelist for additional options)

Mobility



Art. no. 077N
Movable on twin castors
Ø 75 mm, each with total lock



Art. no. 049
Movable on twin castors, Ø 100 mm
with central lock



Art. no. 044N
Wheel-lifting/fixing mechanism
with 4 foot levers for controlled
lowering

Paper roll holder




Art. no. 010F/K
Paper roll holder, mounted under head
part (not protruding), for width 700 mm



Art. no. 010F
Paper roll holder, mounted at foot part
(protruding: 70 mm), for width 700 mm

Nose opening



Art. no. 1615-EF
Loss of nose opening and upholstery
part to close

Add-on elements



Art. no. B1410
1 piece of side rail, length: 350 mm



Art. no. B1430
1 attachment clamp for side rail,
suitable for pipe: Ø 18 mm

Arm rests



Art. no. B1440
1 arm rest for side rail, without
attachment clamp, swivelling and
removable, additional revolving and
tiltable

Side guards



Art. no. 1005
1 side guard to fit on side rail, with
attachment clamp, chromium-plated,
length: 800 mm

Upholstery options

Art. no. B1770
Soft covering for pleasant resting comfort,
for width 700 mm, thickness 80 mm

Art. no. B1772
Upholstery medium-soft (cold foam cover),
for width 700 mm, thickness 70 mm

Updated 04.2023_01
changes concerning technical data are subject to change