


Series 4800E

3-sectioned therapy table with electric motorized height adjustment

Model 4840-00E white powder coated (RAL 9010)
Model 4840-03E white aluminium powder coated (RAL 9006)
Model 4840-04E grey aluminium powder coated (RAL 9007)



German manufacture










All pictures:
Model 4840-00E
with options:
Wheel-lifting/fixing mechanism (045N),
surrounding foot control (095E),
head part 3-sectioned (055-10);
upholstery colour: 64-701 stone blue

Basic equipment:

- Lying surface 3-sectioned
- Multifunctional and independent adjustment of upholstery sections
- Possibilities of positioning the surface parts: Roof shape, Trendelenburg, reverse Trendelenburg
- Table not movable, including height adjustment screws to compensate for uneven ground
- High adjustment range
- Electromotive drive control by foot switch (option: hand switch)
- Integrated safety function for automatic blocking of electromotive drive
- Freewheeling downwards, automatic declutch in case of resistance
- Dewert-hygiene-equipment



Upholstery:

Lying surface 3-sectioned:

- Head part stepless up- (+45°) and downwards (-35°) adjustable by gas spring
- Length head part: 450 mm
- Head part with nose opening and upholstery part to close
- Middle part adjustable to roof position by gas spring
- Length middlepart: 530 mm
- Foot part stepless adjustable (+65°) to sitting position by 2 gas springs
- Length foot part: 1020 mm
- Thickness of upholstery: 60 mm

Technical data:

- Electronical height adjustment: 510 mm to 1030 mm in 18 sec.
  by AC-motor with a pressure capacity of 10000 N
- Width of lying surface: 700 mm
- Length of lying surface: 2000 mm
- Load capacity: 250 kg; according to IEC 60601 tested with fourfold static load

kg

Design:

2 different variants for artificial leather eligible for upholstery covering:


- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



vertical
adjustment

Available as an option (see pricelist for additional options)

Mobility



Art. no. 040N
Movable on comfort twin castors Ø 100 mm, each with total lock, basic height -15 mm



Art. no. 046
Movable on twin castors, Ø 100 mm with central lock, basic height +20 mm



Art. no. 045N
Wheel-lifting/fixing mechanism with 4 foot levers for controlled lowering

DEWERT Präzisions-Stahlrohrmöbel

Additional motor control options



Art. no. 090E
Hand switch (low voltage), with magnet on the back to fasten



Art. no. 091E
Foot switch (low voltage)



Art. no. 092E
1 foot switch fixation mounted at underframe



Art. no. 094E
Foot switch rail, lateral, on both sides, security adjustment:
Press the switch rail down = the table is lifted
Lift the switch rail up = the table is lowered



Art. no. 095E
Foot switch rail, surround, security adjustment:
Press the switch rail down = the table is lifted
Lift the switch rail up = the table is lowered

Paper roll holder



Art. no. 010K
Paper roll holder, mounted under head part (not protruding), for width 700 mm, not in case of 3-sectioned head part (Art. no. 056)




Art. no. 010F
Paper roll holder, mounted at foot part (protruding: 70 mm), for width 700 mm

Head part options




Art. no. 055-10
Head part 3-sectioned, arm rests downwards adjustable from -25° to +45° and removable
(art. no. 010K not possible)
measures:
Head part: L. 500 mm x W. 370 mm
Arm rests: L. 420 mm x W. 160 mm

Arm rests



Art. no. B088
1 arm rest, fixed, swivelling and removable, additional revolving and tiltable

Nose opening



Art. no. 1615-EF
Loss of nose opening and upholstery part to close

Upholstery options

Art. no. 007
Width of upholstery 800 mm

Art. no. B1772
Upholstery medium-soft (cold foam cover), for width 700 mm, thickness 70 mm