


Series 2000XL: Recovery and patient transport tables

2-sectioned recovery and patient transport table with hydrauic height adjustment

Length head part: 550 mm



Model 2220-00XL/H white powder coated (RAL 9010)
Model 2220-03XL/H white aluminium powder coated (RAL 9006)
Model 2220-04XL/H grey aluminium powder coated (RAL 9007)

German manufacture



Model 2220-00XL/H
with option:
Storage basket (1850), picture shows chromium-plated push handles;
upholstery colour
64-717 deep sea

Basic equipment:

- Pedals for foot-controlled height positioning on both sides of the table
- Table movable on 4 comfort twin castors, Ø 125 mm, with central lock, 1 castor with direction lock
- 2 push handles at head and foot part, RAL 9006 powder coated
- 2 side guards, lateral lowerable, chromium-plated, length: 1060 mm
- 4 wall protection bumpers, Ø 200 mm
- Recovery upholstery
- Dewert-hygiene-equipment



Upholstery:

Lying surface 2-sectioned:

- Head part adjustabe from 0° to +30° by 2 metal rachets
- Length head part: 550 mm
- Fixed foot part
- Length foot part: 1400 mm
- Upholstery thickness: 80 mm



Model 2220-00XL/H with option:
Without basic side guards, side guards tiltable (2x1040), paper roll holder (1267K), movable on twin cators with central lock (1180N), picture shows chromium-plated push handles;
upholstery colour: 64-708 cherry

Technical data:

- Manual foot operated height adjustment: 570 mm to 990 mm by hydraulic pump, pressure capacity of 10000 N
- Length of lying surface: 1950 mm | Length overall: 2200 mm
- Width of lying surface: 650 mm | Width overall: 800 mm | optional width: 800 mm, width overall: 950 mm
- Load capacity: 250 kg


kg

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology





Available as an option (see pricelist for additional options)

Mobility



Art. no. 1183
Without basic mobility



Art. no. 1195
Movable on comfort twin castors, Ø 150 mm with central lock, antistatic, basic height + 70 mm



Art. no. 1199
Footbar for central lock, at foot end

without picture

Art. no. 1190
5th castor, connectible, in the centre of the under frame

Side guards



Art. no. 1080-EF
Without basic side guard (1 piece)



Art. no. 1040
1 side guard tiltable, chromium-plated, completely made of steel tube, length: 700 mm, with wall protection bumpers Ø 150 mm



Art. no. 1005
1 side guard to fit on side rail, with attachment clamp, chromium-plated, length: 800 mm

Head part options



Art. no. 1655-30
Comfortable stepless adjustment of head part by gas spring (0° to +60°)

Upholstery options

Art. no. 2220-800
Width of upholstery 800 mm, width overall 950 mm

Art. no. 1765 / 1780
Soft covering for pleasant resting comfort, for width 650 mm / 800 mm, thickness 80 mm

Paper roll holder



Art. no. 1267K / 1268K
Paper roll holder, mounted under head part, for width 650 mm / 800 mm

Add-on elements



Art. no. 1410
1 piece of side rail, length: 350 mm



Art. no. 1430
1 attachment clamp for side rail, suitable for pipe: Ø 18 mm



Art. no. 1450
Infusion pole for side rail, without attachment clamp



Art. no. 1850
Storage basket with mounting, measures:
L. 530 x W. 370 x H. 80 mm

Updated 06.2022_01
changes concerning the technical data are subject to change