
# Series 2000XL: Recovery and patient transport tables

## 2-sectioned Recovery and patient transport table with hydraulic height and inclination adjustment

**Length head part: 750 mm**

**Model 2310-00XL/H white powder coated (RAL 9010)**
**Model 2310-03XL/H white aluminium powder coated (RAL 9006)**
**Model 2310-04XL/H grey aluminium powder coated (RAL 9007)**




German manufacture



Both pictures:
Model 2310-00XL/H
with options: Width of upholstery 800 mm (2310-800), paper roll holder (1268K), side rail (1410), Storage basket (1850), picture shows chromium-plated push handles;
upholstery colour
64-717 deep sea

### Basic equipment:

- Pedals for foot-controlled height positioning on both sides of the table
- Table movable on 4 comfort twin castors, Ø 125 mm, with central lock, 1 castor with direction lock
- 2 push handles at head and foot part, RAL 9006 powder coated
- 2 side guards, lateral lowerable, chromium-plated, length: 1060 mm
- 4 wall protection bumpers, Ø 200 mm
- Recovery upholstery
- Dewert-hygiene-equipment

### Upholstery:

Lying surface 2-sectioned:

- Head part stepless adjustable from 0° to +75° by gas spring (sitting position)
- Length head part: 750 mm
- Fixed foot part
- Length foot part: 1200 mm
- Lying surface incinable till 12° (Trendelenburg) by gas spring
- Upholstery thickness: 80 mm




Picture shows additionally:
Attachment clamp for side rail (1430), infusion pole for side rail (1450) with accessoires

### Technical data:

- Manual foot operated height adjustment: 640 mm to 1060 mm by hydraulic pump, pressure capacity of 10000 N
- Length of lying surface: 1950 mm | Length overall: 2200 mm
- Width of lying surface: 650 mm | Width overall: 800 mm | optional width: 800 mm, width overall: 950 mm
- Load capacity: 250 kg

kg

### Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



vertical
adjustment

## Available as an option (see pricelist for additional options)

DEWERT Präzisions-Stahlrohrmöbel

### Mobility



Art. no. 1183
Without basic mobility



Art. no. 1195
Movable on comfort twin castors, Ø 150 mm with central lock, antistatic, basic height + 70 mm



Art. no. 1199
Footbar for central lock, at foot end

without picture

Art. no. 1190
5th castor, connectible, in the centre of the under frame

### Side guards



Art. no. 1080-EF
Without basic side guard (1 piece)



Art. no. 1040
1 side guard tiltable, chromium-plated, completely made of steel tube, length: 700 mm, with wall protection bumpers Ø 150 mm



Art. no. 1005
1 side guard to fit on side rail, with attachment clamp, chromium-plated, length: 800 mm

### Paper roll holder



Art. no. 1267K / 1268K
Paper roll holder, mounted under head part, for width 650 mm / 800 mm

### Upholstery options

Art. no. 2310-800
Width of upholstery 800 mm, width overall 950 mm

Art. no. 1765 / 1780
Soft covering for pleasant resting comfort, for width 650 mm / 800 mm, thickness 80 mm

### Add-on elements



Art. no. 1410
1 piece of side rail, length: 350 mm



Art. no. 1430
1 attachment clamp for side rail, suitable for pipe: Ø 18 mm



Art. no. 1450
Infusion pole for side rail, without attachment clamp



Art. no. 1850
Storage basket with mounting, measures:
L. 530 x W. 370 x H. 80 mm

Updated 06.2022_01
changes concerning the technical data are subject to change