
Series 2000XL: Recovery and patient transport tables

2-sectioned recovery and patient transport table with hydraulic height and inclination adjustment and detachable mattress


Length head part: 730 mm

Model 2410-00XL/H white powder coated (RAL 9010)
Model 2410-03XL/H white aluminium powder coated (RAL 9006)
Model 2410-04XL/H grey aluminium powder coated (RAL 9007)



German manufacture





All pictures:
Model 2410-00XL/H
with option:
Storage basket (1850),
picture shows
chromium-plated push
handles;
upholstery colour
64-174 sky

Basic equipment:

- Pedals for foot-controlled height positioning on both sides of the table
- Table movable on 4 comfort twin castors, Ø 150 mm, with cental lock, antistatic, 1 castor with direction lock
- Foot bar, chromium-plated, at foot side, facilitates operation of central locking system
- Push handles at head and foot part, RAL 9006 powder coated, at head part: 2 sectioned
- 2 side guards, lateral lowerable, chromium-plated, length: 1060 mm
- 4 wall protection bumpers, Ø 200 mm
- Detachable mattress
- Dewert-hygiene-equipment

Upholstery:

Lying surface 2-sectioned:

- Head part stepless adjustable from 0° to +70° (sitting position) by gas spring
- Length head part: 730 mm
- Fixed foot part
- Length foot part: 1200 mm
- Lying surface incinable till 12° (Trendelenburg) by gas spring
- Special feature: Detachable mattress, see on the right

Technical data:

- Manual foot operated height adjustment: 690 mm to 1110 mm by hydraulic pump, pressure capacity of 10000 N
- Length of lying surface: 1930 mm | Length overall: 2030 mm
- Width of lying surface: 650 mm | Width overall: 800 mm | optional width: 800 mm, width overall: 950 mm
- Load capacity: 250 kg


kg

Design:

2 different variants for artificial leather eligible for upholstery covering:


- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology






Available as an option (see pricelist for additional options)

Mobility



Art. no. 1190
5th castor, connectible, in the centre of the under frame

Paper roll holder



Art. no. 1267K / 1268K
Paper roll holder, mounted under head part, for width 650 mm / 800 mm

Storage basket



Art. no. 1850
Storage basket with mounting, measures:
L. 530 x W. 370 x H. 80 mm

Add-on elements



Art. no. 1410
1 piece of side rail, length: 350 mm



Art. no. 1430
1 attachment clamp for side rail, suitable for pipe: Ø 18 mm



Art. no. 1450
Infusion pole for side rail, without attachment clamp

Upholstery options

Art. no. 2410-800
Width of upholstery 800 mm, width overall 950 mm

Special feature: Mattress:

- Detachable mattress, thickness: 10 cm
- Upholstery made of high-quality resilience (cold) foam, covered with leatherette
- Mattress with velcro fastening and zipper at the bottom
- Support plate for mattress made of HPL panels (resistant to disinfectants)
- Head and foot part hanging together without joints



Updated 06.2022_01
changes concerning the technical data are subject to change