

Recovery and patient transport tables

3-sectioned recovery and patient transport table with fixed height

Length head part: 550 mm

Model 71-00 white powder coated (RAL 9010)
Model 71-03 white aluminium powder coated (RAL 9006)
Model 71-04 grey aluminium powder coated (RAL 9007)

German manufacture



Both pictures:
Model 71-00
with options:
Movable on twin castors with central lock (1180N), adjustment of head part by gas spring (1655), picture shows chromium-plated push handle;
upholstery colour
64-708 cherry


Basic equipment:

- Table movable on 4 castors, Ø 100 mm, each with lock
- Longitudinal brace for additional stability
- 1 push handle, chromium-plated, mounted at foot part (optional at head part possible)
- 2 side guards, lateral lowerable, RAL 9006 powder coated, length: 820 mm
- Recovery upholstery
- Dewert-hygiene-equipment

Upholstery:

Lying surface 3-sectioned:

- Head part upwards (+30°) adjustable by 2 metal rachets
- Length head part: 550 mm
- Fixed middle part
- Length middle part: 450 mm
- Foot part stepless adjustable by gas spring (0° to +32°)
- Length foot part: 950 mm
- Upholstery thickness: 80 mm



Technical data:

- Length of lying surface: 1950 mm | Length overall: 2080 mm
- Width of lying surface: 650 mm | Width overall: 750 mm | optional width: 800 mm, width overall: 900 mm
- Height: 690 mm
- Load capacity: 250 kg kg

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology

Available as an option (see pricelist for additional options)

Mobility



Art. no. 1120-EF
Without basic mobility



Art. no. 1140-11
Movable on 4 castors, Ø 125 mm, each with total lock



Art. no. 1180N
Movable on twin castors, Ø 100 mm with central lock



Art. no. 1185
Movable on twin castors, Ø 125 mm with central lock

Side guards



Art. no. 1088-EF
Without basic side guards



Art. no. 1040
1 side guard tiltable, chromium-plated, completely made of steel tube, length: 700 mm



Art. no. 1005
1 side guard to fit on side rail, with attachment clamp, chromium-plated, length: 800 mm

Head part options



Art. no. 1655-20
Head part stepless adjustable by gas spring (0° to +55°)



Art. no. 1610
Nose opening

Art. no. 1620
Upholstery part to close

Paper roll holder



Art. no. 1210 / 1220
Paper roll holder,
for width 650 mm / 800 mm

Wall protection bumpers

without picture

Art. no. 1172
4 wall protection bumpers Ø 100 mm

Add-on elements



Art. no. 1810
1 additional push handle



Art. no. 1850
Storage basket with mounting, measures:
L. 530 x W. 370 x H. 80 mm



Art. no. 1410
1 piece of side rail, length: 350 mm



Art. no. 1430
1 attachment clamp for side rail,
suitable for pipe: Ø 18 mm



Art. no. 1450
Infusion pole for side rail, without attachment clamp

Upholstery options

Art. no. 71-800
Width of upholstery 800 mm

Art. no. 1765 / 1780
Soft covering for pleasant resting comfort, for width 650 mm / 800 mm, thickness 80 mm

Updated 06.2022_01
changes concerning the technical data are subject to change