

Folding wall tables

2-sectioned folding wall table, turnable along the long side

Length head part: 550 mm

Head part left:
Model 125-00/L white powder coated (RAL 9010)
Model 125-01/L chromium-plated
Model 125-03/L white aluminium powder coated (RAL 9006)
Model 125-04/L silver aluminium powder coated (RAL 9007)

Head part right:
Model 125-00/R white powder coated (RAL 9010)
Model 125-01/R chromium-plated
Model 125-03/R white aluminium powder coated (RAL 9006)
Model 125-04/R silver aluminium powder coated (RAL 9007)




German manufacture

All pictures:
Model 125-01/L

Function:

- Turnable to the wall along the long side upwards or downwards
- Easy operation:
  Swinging out of the feet supported by gas spring
- High stability and less load to the wall:
  In functional state the table stands on its feet



Please note:

- Mounting with standard wall holders is only allowed at a solid wall
- For other wall constructions special wall holders (art. No. 125-1380) with additional stabilization at the floor are necessary
  (possible if there is no baseboard at the table´s location)



Technology:

- Frame construction made of square tube
  30 x 30 x 2 mm
- Galvanised steel hinges connect the upholstery parts
- 2 stable parallel mounted safety metal rachets for positive adjustment of head part

Basic equipment:

- Head part positive adjustable
- Dewert-hygiene-equipment

Included in the scope of supply:
- 2 wall holders for table fixation
- 1 safety holder including belt for safety fixation if the table is turned upwards
- 1 spacer as a swinging stop if the table turned downwards

Upholstery:

Lying surface 2-sectioned:

- Head part upwards (+30°) adjustable by 2 metal rachets
- Length head part: 550 mm
- Fixed foot part
- Length foot part: 1400 mm
- Upholstery thickness: 80 mm

Technical data:

- Length overall: 1950 mm
- Length including wall holders: 2040 mm
- Width overall: 650 mm
- Height: 700 mm
- Weight: 35 kg
- Load capacity (only for solid walls): 200 kg

kg

Design:

2 different variants for artificial leather eligible for upholstery covering:
- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology

Additional option:

without picture

Art. no. 125-1380
Special wall holders

DEWERT Präzisions-Stahlrohrmöbel



Updated 05.2021_01
changes concerning the technical data are subject to change